



ALFAGEN

ALFAGEN LABORATUVAR
MALZEMELERİ İÇ VE DIŞ TİCARET
LTD. ŞTİ.

HÜCRE SAYIM CİHAZI

SON 5 ADET!!!



NanoEnTek

<http://www.nanoentek.com/>

Counting with EVE, affordable automatic cell counting

Key features and benefits

★ **Fine distinction of clumped cell**

An advanced analysis algorithm provides accurate results in clumped cell

★ **Applicable for broad range of cell sizes and types**

From Primary (tissue, blood) cell lines to stem cells

★ **User friendly**

7 inch LCD touch screen Benchtop size No maintenance needed

★ **Data storage and analysis**



Stores 500 test results, transfers the data to PC using USB drive
Provided PC software will give you more in-depth analysis

★ **Cell size gating function**

Select range of cell size



Total cell counting result (Clumped cell)




The MCF7 breast cancer cell line

Clumped cells (HepG2, MCF7, F9) were counted with EVE™, a hemocytometer and the competitor's automated cell counter.

The competitor's cell counter was shown significantly less numbers in clumped cell.

Cell counting results were obtained seven times repetition.

Correlation of EVE™ and manual counting



Measuring from the EVE™ extend further along the high concentration range than hemocytometer readings.

Cell lines validated on EVE™

Cell Type	Animal	Organ	Growth Properties
HeLa	Human	Skin	Adherent
NIH-3T3	Mouse	Embryo	Adherent
U2-OS	Human	Bone	Adherent
Jurkat	Human	Blood	Suspension
KG-1	Human	Blood	Suspension
HepG2	Human	Liver	Adherent
Hep3B	Human	Liver	Adherent
LNCaP	Human	Prostate	Adherent
SH-SY5Y	Human	Brain	Adherent
SCN2.2	Rat	Brain	Adherent
F9	Mouse	Embryo	Adherent
MCF7	Human	Breast	Adherent
A549	Human	Lung	Adherent
GH3	Rat	Pituitary gland	Adherent

In accuracy and precision, EVE™ is comparable to the hemocytometer for all cell lines.

Faster ! Accurate ! Easy to use !

Convenient for use

1. Load the sample



2. Adjust the focus



3. Obtain the results



SON 5 ADET İÇİN **KAÇIRILMAYACAK YENİ YIL** **KAMPANYASI**

SADECE

2.000 \$ + %18 KDV

İzmir Merkez: Tel: 0(232) 375 25 30 , 0(232) 344 25 30 Faks: 0(232) 375 25 31 /
E-Posta: izmir@alfagen.com.tr

Ankara Şube: Tel: 0(312) 286 88 48 Faks: 0(312) 286 88 40
E-posta: ankara@alfagen.com.tr

İstanbul Şube: Tel: 0(533) 778 36 32 / 0(533) 966 51 86
E-posta: istanbul@alfagen.com.tr